



**Agarose, 100g
(QT200-1101)**

Product	Agarose
CAS	9012-36-6
Appearance	White to off-white powder
EEO	≤ 0.13
Gelling Point	36°C±1.5°C (1.5% gel)
Melting Point	88°C±1.5°C (1.5% gel)
Solubility	Clear colorless solution at 1gr in 100ml water
Moisture	≤ 10%
Gel Strength	≥1200 g/cm ² (1% Gel)
Sulfate	≤ 0.15%
Ash	≤ 0.5%
DNase & RNase	None Detected
Protease	None Detected
Endonuclease	None Detected
Storage Temperature	Room Temperature